





Clarity Chromatography Station – Compatibility Table

Compatibility Table

	 Clarity	 Clarity Lite	Note
Windows 7, Vista (both 64 bit)	INT7, INT9, U-PAD2, Net-PAD	INT7, INT9, U-PAD2	<i>From version 3.0</i>
Windows 7 (32 bit)	INT7, INT9*, U-PAD, U-PAD2*, CB20, Net-PAD	INT7, INT9, U-PAD, U-PAD2	<i>From version 2.8</i>
Windows Vista (32 bit)	INT7, INT9*, U-PAD, U-PAD2*, CB20, Net-PAD	INT7, INT9, U-PAD, U-PAD2	<i>From version 2.5</i>
Windows 2000/XP	INT7, INT9*, U-PAD, U-PAD2*, CB20, Net-PAD	INT7, INT9*, U-PAD, U-PAD2*	<i>From version 2.4.4</i>
Windows 98/ME	INT7, U-PAD, INT5, CB20, Net-PAD	INT7, U-PAD	<i>Supported only until version 2.4.1</i>

* INT9, UPAD2 from version 2.4.4

- HW dongle required for all systems from version 2.4.0

PC Configurations

	Minimal	Recommended
Windows 7	see Windows 7 specification	see Windows 7 specification
Windows Vista	see Windows Vista specification	see Windows Vista specification
Windows XP/2000	PC Pentium III/700 MHz, 256 MB RAM	PC Pentium 4/2 GHz, 512 MB RAM
Windows 98/ME	PC Pentium/200 MHz, 32 MB RAM	PC Pentium/450 MHz, 128 MB RAM
Monitor	Resolution 1024x768, 64K (16 bit High color)	Resolution 1280x1024, 64K (16 bit High color)

- USB or LPT for a dongle
- Free PCI (ISA) slot for internal A/D card or USB port for UPAD or LAN port for Net-PAD
- CD ROM drive